



Nala™ patient chair





Soft front edge on seat and soft, extended arm pads help patients guide and stabilize themselves during transfer.



Barrier-free foot area gives patients space to position their feet for leverage during egress.



Nala Harmonic™ Tilt and Ottoman
CR43/CR900
W 33 H 50.5 D 28.5 / W 25.25 H 16.25 D 16

Nala


patient chair

Science of seating, comfort of sitting.

Patient comfort takes form in the Nala chair. Its design and construction incorporate ergonomic principles and technologies from seating leader Herman Miller and its award-winning chair lineup. Its function and materials recognize the realities of patient care environments. The result: a chair that works with the patient's body. Nala's tilt, support, and contours all follow the movement and form of the human body.




Enhanced Features



Ergonomic Design

Nala's Harmonic tilt mimics the natural movement of the patient's body. The chair's back and seat move synchronously, opening up as the patient reclines and tilting around the body's natural pivot points. Patients can rest anywhere within a 24-degree recline range. Dampening cylinder controls movement for heavy or light patients.



Support + Comfort

Nala's seat-and-back suspension is a series of firm, flexible straps encased in a layer of high-density foam. The combination offers a perfect balance: responsive and even support over the entire chair surface and a soft outer layer for a comfortable sit and feel.




Environmental Design

Nala achieved McDonough Braungart Design Chemistry (MBDC) Silver Cradle to Cradle™ (C2C) certification, the first patient chair to receive this certification. It follows environmental protocols established by Herman Miller, in conjunction with MBDC and its C2C design protocols. Nala has 33 percent recycled content, and is 86 percent recyclable.

New Level of Performance



Nala Harmonic Tilt (arm raised)
CR43
W 33 H 50.5 D 28.5



Pivoting arms give caregivers greater access to patients during transfers and treatments.

Materials Options



Nala's upholstery options are designed to meet the performance and durability standards of healthcare organizations. Environmentally friendly, proprietary textile options are also available, including PVC-free polyurethanes and antimony-free polyesters.

Nala offers designers choices for combining textiles and colors, including the five upholstery configuration options shown here.

Upholstery configurations



Configuration 1



Configuration 2



Configuration 3



Configuration 4



Configuration 5

Design Features



Recline and arm pivot controls on each arm give caregivers and patients control from either side of the chair and with either hand.



Leaf spring in tilt mechanism acts as a power source that minimizes the physical work of moving seated patients forward.



Integral handle assists caregivers in moving chair.



Wheels on rear legs engage when chair is tipped back.



Headrest pillow adjusts 9 inches. Edges curve to support the patient's head.

Statement of Line



Nala
Harmonic Tilt



Nala
Flex Back



Nala
Fixed Back



Nala
Ottoman

Notes

- Featured textiles: Herman Miller Quilty, Squash; Maharam Droplet, Cayenne; Maharam Chime, Flame; Maharam Kernel, Shilling
- Nala supports users who weigh up to 350 pounds
- 12-year, 24/7 warranty, including parts and labor
- Designed by Continuum, Boston, Massachusetts
- Contact Brandrud Customer Service at 253.838.6500 for detailed cleaning guidelines and suggested cleaning solutions