

sittris

Sittris PA



Introducing Sittris PA

An effective patient chair must be much more than simply a comfortable place to sit. It needs to address the issues faced daily by healthcare workers and their patients and provide attributes that encourage patient recovery:

- providing gentle, relaxing therapeutic movement to enhance patient convalescence
- providing pressure relief in the seat area to stimulate blood circulation – reducing the risk of skin breakdown, particularly in areas of bony protrusions such as the coccyx
- providing comfortable support for patients' arms without putting unnecessary pressure on the elbows – reducing the chance for elbow skin shear
- providing nursing staff easy access to patients for chair-side treatments and allow for safe, discreet placement of IV poles, oxygen equipment, etc.
- providing comfortable leg positioning while reducing downward pressure on patients' heels – reducing the risk of skin breakdown
- providing easy to clean surfaces that resist bacterial, viral and microbial growth

Of course, an effective patient chair must be comfortable as well – offering superb lumbar support and effectively “nestling” all patients regardless of height allowing them to relax and “snooze.”

We designed Sittris PA with all this in mind – a chair truly designed for patient care.

Designed for patient care

Back rest

Both the knit and fully upholstered high back design offer superb lumbar support and effectively “nestle” all patients regardless of height.

Head rest

The optional headrest supports the neck and cradles the patient's head further helping them to relax and “snooze”.



Arm design

Both the urethane and fully upholstered arms are designed to provide comfortable support for the patient's arms without putting unnecessary pressure on the elbows, reducing the chance for elbow skin shear.



Flex-i-back

The flex-i-back version provides gentle, relaxing therapeutic movement to enhance patient recovery.

Patient Access

The Access Design Back provides nursing staff easy access to patients for chair side treatments and allows for safe, discreet placement of IV poles, oxygen equipment, etc.

Clean Seating

Like all Sittris products, the Sittris PA is designed to allow for easy cleaning and all Sittris textiles are cleanable with soap and water. In addition, Sittris textiles are anti-bacterial, anti-microbial, have a moisture barrier, but are also breathable for patient comfort.



Integrated pressure relief

The pressure relief system that is built into the seat stimulates blood circulation to reduce the risk of skin breakdown, particularly in areas of bony protrusions such as the coccyx.

Small footprint

The sleek, thin design and small footprint of the Sittris PA is designed to maximize comfort, while offering efficient space utilization.

Ottoman

The Sittris Ottoman is designed for comfortable leg positioning while reducing downward pressure on the patient's heels, reducing the risk of skin breakdown.





SITRIS PA KNIT BACK with Urethane Arms



SITTRIS PA with Upholstered Back and Arms

Sittris PA

Designed by Kerr and Company

Sittris PA Knits



Sunset



Autumn



Sahara



Osoyoos



Savannah



Shenandoah



Alleghany



Gobi



Okanagan



Adriatic



Fjord



Everglades

sittris

Sittris
1406 Birchmount Road
Toronto, Ontario
Canada M1P 2E3
www.sittris.com

416 701 4673
877 701 4673